




# HUMAN BRAIN ORGANOID-ON-CHIP PLATFORM TO IMPROVE ORGANOID REPRODUCIBILITY AND SCALABILITY FOR PHARMACEUTICAL STUDIES




Héloïse Castiglione<sup>1,2,3</sup>, Lucie Madrange<sup>2,3</sup>, Camille Baquerre<sup>1</sup>, Benoît G. C. Maisonneuve<sup>1</sup>, Johan Renault<sup>1</sup>, Frank Yates<sup>2,3</sup>, Pierre-Antoine Vigneron<sup>2,3</sup>, Jessica Rontard<sup>1</sup>, Thibault Honegger<sup>1</sup>  
<sup>1</sup> NETRI, 69007 Lyon, France, <sup>2</sup> SupBiotech, Ecole d'ingénieurs en Biotechnologies, 94800 Villejuif, France, <sup>3</sup> Service d'Etude des Prions et des Infections Atypiques (SEPIA), Institut François Jacob, Commissariat à l'Energie Atomique et aux Energies Alternatives (CEA), Université Paris Saclay, 92260 Fontenay-aux-Roses, France

## BACKGROUND

### Cerebral organoid



### Microfluidic device (NETRI)



Cerebral organoids are emerging as relevant alternatives for modeling human brain cellular organization and development in 3D. To facilitate the use of organoids for large-scale drug screening, they require a gain of reproducibility and scalability (Castiglione *et al.*, 2022). To address this challenge, we have:

- Developed a Brain Organoid-on-a-Chip platform, combining cerebral organoid culture within a NETRI microfluidic device
- Optimized an on-chip culture protocol for cerebral organoids, enhancing their reproducibility
- Established a scoring system & a prediction algorithm, for quality control and toxicological evaluation of cortical organoids exposed to chemical compounds

Using this platform, we have initiated neurotoxicity evaluations of vanillin and biphenyl-2-ylamine.

## EXPERIMENTAL DESIGN

### DUPLEX WELL MICROFLUIDIC DEVICE


#### Adapted to 3D cell culture:

Two compartments separated by a porous membrane:

- Open well for 3D culture
- Perfusion channel

#### Adapted to industrial transfer

#### Pumpless




Duplex Well (©NETRI)

### BRAIN ORGANOID-ON-CHIP PROTOCOL & NEUROTOXICITY STUDIES

#### Standardized on-chip protocol for cortical organoid culture

- Viability and expected morphology up to D+120
- Expected RNA expression levels at D+60
- Improved cytoarchitectural organization at D+60
- Anti-adherence protocol to limit organoid adherence onto membrane



On-chip culture conditions and neurotoxicity evaluations.  
Cortical organoid generation and culture protocol adapted from Xiang *et al.*, 2019; hiPSCs: human induced pluripotent stem cells, EB: embryo body.


Vanillin (CAS n°121-33-5): 100, 1 000, 10 000 nM  
Biphenyl-2-ylamine (CAS n°90-41-5): 20, 200, 2 000 µM (meOH)

## RESULTS


### IMPROVED REPRODUCIBILITY ON-CHIP COMPARED WITH CONVENTIONAL SUPPORT

#### More reproducible sizes of cerebral organoids on-chip compared to controls cultured in 24-well plates

- Lower size dispersions for organoids on-chip compared with organoids cultured in 24-well plates → reduced variability
- Smaller average organoid sizes on-chip could be due to the culture in a constraint environment
- Observed with two hiPSCs lines



Cerebral organoid sizes and growth evolutions from D+18 until D+60 of culture, for on-chip organoids and control organoids cultured in 24-well plates. On-chip organoids exhibit more reproducible surface areas (24-well plate: n=9; Duplex Well: n=27).



Immunofluorescence staining of neural progenitors (SOX2) and neurons (TUBB3) in cerebral organoids after 60 days of culture. (Thunder microscope, Leica, objective 10X).


#### More reproducible cytoarchitectural organization in organoids on-chip than in 24-well plates

Rosettes with expected morphologies (oval shapes) & well-organized patterns (present along organoid borders)

### NEUROTOXICITY STUDIES USING THE BRAIN ORGANOID-ON-CHIP PLATFORM


#### Example 1: acute exposure (10 000 nM vanillin)

##### Optimal morphology




Morphology of cortical organoids at day 61: examples of an unexposed control organoid and vanillin-exposed organoid #1 (brightfield, 5X).

##### Similar growth profile compared to controls



Cortical organoids growth curves and slopes between exposure (D+41) and D+61 (controls: mean ±SEM, n=4).

##### Expected cell types and optimal cytoarchitectural organization




Immunofluorescence staining of neural progenitors (SOX2), neurons (TUBB3, MAP2), and astrocytes (GFAP) (Thunder microscope, Leica, 20X).

#### General conclusions using the scoring system & prediction algorithm

Vanillin exposures: no discernable impact on morphology, cytoarchitectures & viability


#### Example 2: acute exposure (2 000 µM biphenyl-2-ylamine)

##### Altered morphology




Morphology of cortical organoid #2 at 61 days of culture (brightfield, 5X, circle: cyst).

##### Higher apoptosis & DNA damage levels



Immunofluorescence staining of apoptosis (CCASP3) and DNA damage (γH2AX) markers (Thunder microscope, Leica, 20X).

##### Expected cell types with disorganized cytoarchitectures



Immunofluorescence staining of neural progenitors (SOX2), neurons (TUBB3, MAP2), and astrocytes (GFAP) (Thunder microscope, Leica, 20X).


- Altered pattern of rosettes
- Presence of a large necrotic core

Biphenyl-2-ylamine exposures: altered morphology & disorganized cytoarchitectures in a dose-response manner

## CONCLUSION & PERSPECTIVES

- Optimized cerebral organoid culture protocol on-chip
- Improved intra- and inter-batch reproducibility: in terms of size, growth profile and cytoarchitectural organization
- Brain Organoid-on-Chip platform + scoring system + prediction algorithm: adapted to neurotoxicity evaluations


- Evaluation of the neurotoxicity of other compounds to refine the prediction algorithm
- Modelling the blood-brain barrier
- Our Brain Organoid-on-Chip platform paves the way for drug screening and toxicological assessments



## CONTACT

www.netri.com  
Phone: +33 4 78 23 08 66  
Email: contact@netri.com

www.recherche.supbiotech.fr  
Email: pierre-antoine.vigneron@supbiotech.fr



Presented on 2024-06-10 - 14  
[Seattle, Washington, USA]