

EzWay™ One-Step RT-PCR MasterMix

1. Catalog No.

Cat. No	Product	Size
K0569600	EzWay™ One-Step RT-PCR MasterMix	1ml

2. Storage

1 year at -20°C
 (The product is able to be shipped on blue ice and should be stored immediately at -20°C.)

3. Contents

Component	Cat. No	Packing	Comments
2X EzWay™ One-Step RT-PCR MasterMix	K0569600	1.0ml	Mixture of Hot Taq DNA polymerase, Reverse transcriptase, optimized PCR reaction buffer (MgCl ₂), dNTP, stabilizer, dye and additive
4X Magic Buffer	K0561031	1.0ml	<u>Only use for High GC content</u>

4. Description

EzWay™ One-Step RT-PCR MasterMix is designed for sensitive and specific amplification in one-tube reaction from total RNA transcripts.

EzWay™ One-Step RT-PCR MasterMix is a premixed buffer solution containing MMLV Reverse Transcriptase and Hot Taq DNA Polymerase, dNTPs, tracking & loading dye and stabilizers. The solution contains the proper concentration of MgCl₂ and dNTPs for the better specificity and sensitivity and supplied as 2X concentrated type for convenience.

Simply add primers and RNA template diluted in DEPC-treated water into equal volume of EzWay™ One-Step RT-PCR MasterMix. Pipetting steps are minimized, reducing the possibility of errors and contamination

- Synthesis of cDNA and sequentially PCR amplification from target gene in one-tube.
- Specialized buffer composition for highly sensitive and specific RNA amplification
- Simple set-up and minimal optimization required

5. Application

- Rapid detection of target gene transcript from total RNA or purified mRNA

6. PCR Amplification

This protocol serves as a guideline for one-step RT-PCR. Reverse transcription and PCR are carried out sequentially in the same tube. Optimal reaction conditions, such as incubation times, temperatures, and amount of template RNA, may vary and must be individually determined.

Important notes before starting

- EzWay™ One-Step RT-PCR MasterMix requires initial activation by incubation at 95°C for 15 min before amplification can take place. This incubation also inactivates the reverse transcriptases. Do not heat activate Hot Taq DNA Polymerase until the reverse transcriptase reaction is finished.
- EzWay™ One-Step RT-PCR MasterMix is designed to be used with gene-specific primers at a final concentration of 0.2-0.6 µM. The use of random oligomers or oligo-dT primers is not recommended since it will result in the amplification of nonspecific products.
- Set up all reactions on ice.
- Make sure the thermal cycler is preheated to 42-55°C before placing samples in it.
- EzWay™ One-Step RT-PCR MasterMix contains the proper conc. of MgCl₂ which will produce satisfactory results in most cases. However, if a higher Mg²⁺ concentration is required, increase MgCl₂ conc. with the higher conc. of MgCl₂ solution.
- RNase-free environment should be maintained during RNA isolation and reaction setup.
- Set up the reaction mixtures in an area separate from that used for RNA preparation or PCR product analysis.
- Use disposable tips containing hydrophobic filters to minimize cross-contamination.

1. Thaw primer solutions

Keep on ice after complete thawing, and mix well before use.

Optional: Prepare a primer mix of an appropriate concentration (see below) using the RNase-free water. This is recommended if several amplification reactions using the same primer pair are to be performed. The final volume of diluted primer mix plus the template RNA, added at step 4, should be 25 µl per reaction.

Table 1.

Component	Final Concentration	Volume/reaction
2X EzWay™ One-Step RT-PCR Universal qMasterMix	1X	25 µL
5' Primer	0.2-0.6 uM	Variable
3' Primer	0.2-0.6 uM	Variable
Distilled water	-	Variable
Template RNA	total RNA : 1 ng ~ 2 µg Poly(A)+ RNA : 10 pg ~ 500 ng	Variable
Optional) Magic Buffer	5-12.5ul/reaction	Variable
Optional) RNase inhibitor	5-10 units/reaction	Variable
Total reaction volume		50 µL

- Mix by vortexing briefly, and dispense 25 µl into each PCR tube according to Table 1.**

It is important to mix EzWay™ One-Step RT-PCR MasterMix before use to avoid localized differences in salt concentration. EzWay™ One-Step RT-PCR MasterMix is provided as a 2X concentrate (i.e., a 25-µl volume of the EzWay™ One-Step RT-PCR MasterMix is required for amplification reactions with a final volume of 50 µl). For volumes smaller than 50 µl, the 1/1 ratio of EzWay™ One-Step RT-PCR MasterMix to diluted primer mix and template should be maintained as defined in Table 1. A negative control (without template RNA) should be included in every experiment. It is recommended that the PCR tubes are kept on ice until they are placed in the thermal cycler.

- Distribute the appropriate volume of diluted primer mix into the PCR tubes containing EzWay™ One-Step RT-PCR MasterMix.**
- Add template RNA to the individual PCR tubes.**
- Program the thermal cycler according to the program outlined in the table below.**

Step		Temp.	Time	Cycles
Reverse Transcription		42-55°C	30-60 min	1
Initial Denaturation		95°C	15 min	1
Cycling	Denaturation	94°C	0.5–1 min	25-40
	Annealing	50-68°C	0.5–1 min	
	Extension	72°C	1 min	
Final extension		72°C	10 min	1

Note: Pre-heating is required to make DNA polymerase active. Reverse transcriptase is inactivated and the cDNA template is denatured.

- Start the RT-PCR program while PCR tubes are still on ice. Wait until the thermal cycler has reached 50°C. Then place the PCR tubes in the thermal cycler.**

Note: After amplification, samples can be stored overnight at 2–8°C, or at –20°C for longer-term storage.