

Goat Poly Anti-Mouse IgG (H&L), HRP

- 1. Catalog No.** K0211589
- 2. Quantity** 1 mg/ml, 1 ml
- 3. Storage & Stability** Store at 4°C
- 4. Description** Goats were immunized with highly purified Mouse IgG isolated from pooled normal mouse serum. After multiple immunizations in Freund's adjuvant the Goats were bled and serum collected. The anti IgG(H&L) chains were immunoaffinity purified off an IgG-containing solid phase immunoabsorbant and conjugated with horseradish peroxidase (HRP).
- 5. Specificity** This antibody reacts with mouse IgG and with light chains common with other mouse immunoglobulins as determined by IEP and ELISA techniques. It is suitable for blotting, ELISA and IHC applications. This antibody may cross react with IgG from other species.
- 6. Source** Goat
- 7. Immunogen** Mouse IgG (H&L)
- 8. Formulation** This antibody is provided in buffer containing HEPES (-Saline with Stabilizer) with 0.002% Thimerosal
- 9. Purity** N/A
- 10. Application**
- ELISA Detection**
Working Dilution : 0.2 ug/ml
- Western Blotting**
Working Dilution : 0.5 ug/ml
- Investigator should determine optimal working dilutions experimentally. Prepare working dilution immediately before use.