

---

## Goat Poly Anti-Rabbit IgG (H&L), HRP

- 1. Catalog No.** K0211708
- 2. Quantity** 1.0 mg/ml, 1 ml
- 3. Storage & Stability** Store at 4°C
- 4. Description** Goats were immunized with purified Rabbit IgG isolated from pooled normal rabbit serum. After multiple immunizations in Freund's adjuvant the Goats were bled and serum collected. The anti IgG(H&L) chains were immunoaffinity purified off an IgG-containing solid phase immunoabsorbant and conjugated with horseradish peroxidase (HRP).
- 5. Specificity** This antibody reacts with rabbit IgG and with light chains common with other rabbit immunoglobulins as determined by IEP and ELISA techniques. It is suitable for blotting, ELISA and IHC applications. This antibody may cross react with IgG from other species.
- 6. Source** Goat
- 7. Immunogen** Rabbit IgG
- 8. Formulation** This antibody is provided in buffer containing HEPES (-Saline with Stabilizer) with 0.002% Thimerosal.
- 9. Purity**
- 10. Application** **ELISA detection** : 0.05 ug/ml  
**Western Blotting** : 0.2 ug/ml
- Investigator should determine optimal working dilutions experimentally. Prepare working dilution immediately before use.