

## ELISA well plate Sealer

**1. Catalog No.** K0331010

**2. Quantity** 100 EA

**3. Storage** Store at R.T.

**4. Description** Designed for in vitro diagnostic assays

- Standard adhesive-backed microplate sealer suitable for multi well plate sealing where protection of sample contamination and evaporation prevention are desired.
- Specification

Matreial	PET
Sterility	Non-sterile
Size	13 x 8 cm
Thickness	0.05 mm
Temperature Tolerance range	0-40°C
Color	Clear

- Recommended use for 96 well micro plate.

**5. Contents**

Description	Amount
ELISA well plate sealer	100 EA

**6. Application** ELISA