

## EzWay™ Coating Buffer (pH 9.6)

**1. Catalog No.** K0331022

**2. Quantity** 500 ml

**3. Storage** Store at 2-8°C. Do not freeze.

**4. Description** ■ Coating buffer for ELISA assay.

**5. Contents**

Description	Amount
Coating Buffer (pH 9.6)	500 ml

**6. Coating**

- Coating Solution: Resolve the coating material (antigen or antibody) in the coating buffer to make 1 ug/ml (1-10 ug/ml).
- Dispense 100 ul (50-200 ul) of prepared Coating Solution to each well.
- Incubate for overnight at 4°C.