
Blocking Solution (0.1% Casein/PBS)

- 1. Catalog No.** K0331031
- 2. Quantity** 125 ml x 2
- 3. Storage** Store at 2-8°C. Do not freeze.
- 4. Description** This Blocking solution can used in ELISA to blocking plate and Antigen, Antibody, samples. As a blocking solution it binds to any remaining protein binding sites after the primary antigen or antibody is bound to the solid phase surface. As a diluent for antigens, antibodies, and other proteins. it reduce non-specific interaction of ay reactants with each other and with the solid phase surface.
- 5. Contents**
- | Description | Amount |
|-------------------------------------|------------|
| Blocking Solution (0.1% Casein/PBS) | 125 ml x 2 |
- 6. Application** ELISA :. Sample dilution or Antigens, Antibodies dilution.