

EzWay™ pink-ONE TMB Substrate Solution

- 1. Catalog No.** K0331070
- 2. Quantity** 100 ml (1000 well)
- 3. Storage** Store at 2-8°C. Do not freeze.
- 4. Description**
- One-step prestained TMB
 - No mixing of TMB substrate and substrate solution required
 - Pink-color coded indicator : Convenient to track for reagent supply
 - No interference of color to activity or sensitivity of reaction (H_2O_2 makes color disappeared)

5. Kit Contents

| Description | Amount |
|---------------------------------|--------|
| Pink-ONE TMB Substrate Solution | 100 ml |

6. Color Reaction and Reading

- Dispense 100 ul Pink-ONE TMB Substrate Solution into each well.
- After sufficient color development (5-30minutes at room temperature or at 37°C), add 100 ul Stop Solution (2M H_2SO_4 or 0.2M H_2SO_4) to each well.
- Read plates in a microwell plate reader at wavelength setting of 450nm