


# Cytokine ELISA Combo kit

- 1. **Catalog No.** K033KIT-02
- 2. **Quantity** 96 Tests (48 wells x 2 different cytokines)
- 3. **Storage** Store at 4°C.
- 4. **Description**



## 5. Kit Contents

Common Component	Description		Amount
	Choice 1	Choice 2	
Pre-Coated 96 well ELISA microplate	Mouse IFN-gamma Antibody coated	Mouse IL-12 Antibody coated	48 wells x 2
Detection Antibody	Biotinylated Mouse IFN-gamma Antibody	Biotinylated Mouse IL-12 Antibody	1 EA each
Standard Protein	Mouse IFN-gamma Standard	Mouse IL-12 Standard	1 EA each
Color Development Enzyme	Streptavidin-HRP conjugate (600 ul)		1 EA
pink-ONE Assay Diluent	1% BSA in PBS (50 ml)		1 EA
Prestained Color Development Reagent	TMB solution (10 ml)		1 EA
Stop Solution	2M H <sub>2</sub> SO <sub>4</sub> (10 ml)		1 EA
Wash Buffer (20X)	Concentrated PBST for 1L (50ml)		1 EA
Plate Sealer			3 EA

- 6. **ELISA Protocol** Refer to the each cytokine kit manual attached.