


## Cytokine ELISA Combo kit

- 1. **Catalog No.** K033KIT-02
- 2. **Quantity** 96 Tests (48 wells x 2 different cytokines)
- 3. **Storage** Store at 4°C.
- 4. **Description**



### 5. Kit Contents

Common Component	Description		Amount
	Choice 1	Choice 2	
Pre-Coated 96 well ELISA microplate	Cytokine 1 Antibody coated	Cytokine 2 Antibody coated	48 wells x 2
Detection Antibody	Biotinylated Cytokine 1 Antibody	Biotinylated Cytokine 2 Antibody	1 EA each
Standard Protein	Cytokine 1 Standard	Cytokine 2 Standard	1 EA each
Color Development Enzyme	Streptavidin-HRP conjugate (600 ul)		1 EA
Assay Diluent	1% BSA in PBS (50 ml)		1 EA
Prestained Color Development Reagent	TMB solution (10 ml)		1 EA
Stop Solution	2M H <sub>2</sub> SO <sub>4</sub> (10 ml)		1 EA
Wash Buffer (20X)	Concentrated PBST for 1L (50ml)		1 EA
Plate Sealer			3 EA

- 6. **ELISA Protocol** Refer to the each cytokine kit manual.