


# Cytokine ELISA Combo kit, pink-ONE

- 1. **Catalog No.** K033KIT-02P
- 2. **Quantity** 96 Tests (48 wells x 2 different cytokines)
- 3. **Storage** Store at 4°C.
- 4. **Description**



## 5. Kit Contents

| Common Component                     | Description                               |                                   | Amount       |
|--------------------------------------|-------------------------------------------|-----------------------------------|--------------|
|                                      | Choice 1                                  | Choice 2                          |              |
| Pre-Coated 96 well ELISA microplate  | Mouse MIP-1 beta Antibody coated          | Mouse IP-10 Antibody coated       | 48 wells x 2 |
| Detection Antibody                   | Biotinylated Mouse MIP-1 beta Antibody    | Biotinylated Mouse IP-10 Antibody | 1 EA each    |
| Standard Protein                     | Mouse MIP-1 beta Standard                 | Mouse IP-10 Standard              | 1 EA each    |
| Color Development Enzyme             | Streptavidin-HRP conjugate (600 ul)       |                                   | 1 EA         |
| pink-ONE Assay Diluent               | 1% BSA in PBS (50 ml)                     |                                   | 1 EA         |
| Prestained Color Development Reagent | pink-ONE TMB solution (10 ml)             |                                   | 1 EA         |
| Stop Solution                        | 2M H <sub>2</sub> SO <sub>4</sub> (10 ml) |                                   | 1 EA         |
| Wash Buffer (20X)                    | Concentrated PBST for 1L (50ml)           |                                   | 1 EA         |
| Plate Sealer                         |                                           |                                   | 3 EA         |

- 6. **ELISA Protocol** Refer to the each cytokine kit manual attached.