


# Immunoglobulin ELISA Combo kit, pink-ONE

- 1. **Catalog No.** K323KIT-02P
- 2. **Quantity** 96 Tests (48 wells x 2 different immunoglobulins)
- 3. **Storage** Store at 4°C.

4. **Description**



5. **Kit Contents**

Common Component	Description		Amount
	Choice 1	Choice 2	
Pre-Coated 96 well ELISA microplate	Mouse IgE Antibody coated	Mouse IgG Antibody coated	48 wells x 2
Detection Antibody	HRP-conjugated Mouse IgE Antibody	HRP-conjugated Mouse IgG Antibody	1 Vial each
Standard Protein	Mouse IgE Standard	Mouse IgG Standard	1 Vial each
Assay Diluent	0.1% Casein in PBS (30 ml)	1% BSA in PBS (30 ml)	1 EA
Prestained Color Development Reagent	pink-ONE TMB solution (10 ml)		1 EA
Stop Solution	2M H <sub>2</sub> SO <sub>4</sub> (10 ml)		1 EA
Wash Buffer (20X)	Concentrated PBST for 1L (50ml)		1 EA
Plate Sealer			4 EA

- 6. **ELISA Protocol** Refer to the each cytokine kit manual attached.