


Immunoglobulin ELISA Combo kit, pink-ONE

1. **Catalog No.** K323KIT-02P
2. **Quantity** 96 Tests (48 wells x 2 different immunoglobulins)
3. **Storage** Store at 4°C.

4. Description



5. Kit Contents

Common Component	Description		Amount
	Choice 1	Choice 2	
Pre-Coated 96 well ELISA microplate	Immunoglobulin 1 Antibody coated	Immunoglobulin 2 Antibody coated	48 wells x 2
Detection Antibody	HRP-conjugated Immunoglobulin 1 Antibody	HRP-conjugated Immunoglobulin 2 Antibody	1 Vial each
Standard Protein	Immunoglobulin 1 Standard	Immunoglobulin 2 Standard	1 Vial each
Assay Diluent	0.1% Casein in PBS (30 ml)	1% BSA in PBS (30 ml)	1 EA
Prestained Color Development Reagent	pink-ONE TMB solution (10 ml)		1 EA
Stop Solution	2M H ₂ SO ₄ (10 ml)		1 EA
Wash Buffer (20X)	Concentrated PBST for 1L (50ml)		1 EA
Plate Sealer			4 EA

6. **ELISA Protocol** Refer to the each cytokine kit manual.